

# How to use an ELSD



Model 200S

Model 100



INTRODUCING  
Model 2300  
High Sample Throughput.  
Sleek Modern Design.



Model 300S



Model 1300


# Model 2300 Operation



Controlled via six touch sensitive keys  
Hidden behind the molded front

Power  
Home  
<= Left  
=> Right  
Enter  
Next

## Tiered Menu Structure



# Model 1300 Operation

- Controlled via six touch sensitive keys

PWR (power),

A/Z

ESC(escape),

▲ (up),

▼ (down)

ENT (enter)

- Three Tiered Menu



# Model 300S Operation

Only 4 Keys  
for Detector Control

Three Tiered Menu




# Model 200S Operation

Only 2 Keys for Detector Control  
Power Key is used to enter menu  
Autozero Key is used to adjust  
parameters



# PC Control Program



USB Connection to ELSD

# Success with an ELSD

- Clean Mobile Phase Solvents & Modifiers
- Clean gas
- Free flow of exhaust
- Optimized detector temperatures
- Low mobile phase flow rates




# Volatile Mobile Phases


Solvent	Type	% Residue on Evaporation
2-propanol	Reagent Grade	<0.0001
	Spectrophotometric grade	0.0002
Acetonitrile	ACS Grade	<0.001
	HPLC Grade	0.0001
Methanol	Reagent Grade	<0.1
	HPLC Grade	<0.0001
Water	Reagent Grade	<0.1
	HPLC Grade	<0.0001

# Solvent Quality Effects Baseline Noise

HPLC Grade MEOH




Reagent Grade MEOH




# Water Quality Effect on Baseline Noise

HPLC Grade



Reagent Grade



# Volatile Buffers

LCMS Grade is best

	Typical Concentration	pH	Buffering Range
Formic Acid	0.1%	2.7	
Acetic Acid	0.1%	3.3	
Trifluoroacetic Acid	0.1%		
Ammonium Formate	5 – 10mM		2.7 – 4.7
Ammonium Acetate	5 – 10mM		3.7 – 5.7
Ammonium Carbonate	5 – 10mM		6.6 – 8.6

# Starting up an ELSD

- Unpack the instrument and find a spot on the bench with adequate ventilation
- Plug it in but DO NOT TURN IT ON yet.
- If using the PC Control program, Start Windows on the control computer
- Start the control program for the ELSD
- Turn the ELSD on, and select the Universal method from the front panel.
- Connect the liquid inlet to appropriate spot on the HPLC system.
- Fill the drain with mobile phase
- Connect the exhaust to a vent
- Pump mobile phase for a few minutes to wash out the instrument
- Inject your first peak

# ELSD Set Points

## Models 1300 and 300S

### **Universal Conditions**

Spray Chamber 30C  
Drift Tube Temperature 60  
Filter BLT 5

### **Aqueous Conditions**

Spray Chamber 10C  
Drift Tube Temperature 50  
Filter BLT 5

### **Organic Conditions**

Spray Chamber 60C  
Drift Tube Temperature 70  
Filter BLT 5

## Model 2300

### **Universal Conditions**

Spray Chamber 10C  
Drift Tube Temperature 65  
Filter FLT 2

### **Hi FLO Aqueous**

Spray Chamber 20C  
Drift Tube Temperature 85  
Filter FLT 2

### **Low Flo Organic**

Spray Chamber 55C  
Drift Tube Temperature 70  
Filter FLT 5, Low Gain

### **Volatiles (with He)**

Spray Chamber 10C  
Drift Tube Temperature 45  
Filter BFT1

## **When Optimization may be necessary**

- The mobile phase may be difficult to evaporate, high aqueous or volatile buffer content or high flow rates
- The required detection limit is low, 100ng or lower
- The analyte is semi-volatile.

# Parameters to Adjust

## Optimize Conditions for Best Signal to Noise

### Spray Chamber Temperature

Sub ambient for hard to evaporate MP  
Elevated for highly organic MP

### Drift Tube Temperature

Hot enough to evaporate MP  
Not too hot to volatilize analyte

### Filter

Off (or 0) to 10  
Baseline Filter (BLT) vs. RC Filter (FLT)


### Gain

Normal vs EDR (Models 1300, 300S)  
Normal vs. Low (Model 200S and 100, 2300)

### Full Scale

2300 and 1300: 1V or 10mV,  
300S, 200S, 100: 5V or 10mV  
Data systems that require 1V: Change Calibration to 20% at 5V


## Drift Tube Temperature Optimization



In this example a difficult to evaporate mobile phase was selected, 1ml/min Water. The Spray Chamber was set at 30 and the drift tube was varied. At time 0 the drift tube was 60C. The drift tube cooled to 50C by time 5 and returned to 60 after 6 minutes.




## Spray Chamber Temperature Optimization




In this example an easy to evaporate mobile phase was selected, 0.5ml/min Methanol. The Drift Tube temperature was held constant at 40C. The Spray chamber temperature was increased from 10C to 60C. You can see the effects of reducing the spray chamber temperature in the image below.

## Semi-Volatile Analyte Detection




This example show that a small increase in drift tube temperature will decrease the signal of a semi-volatile analyte significantly.

## Filter Type



In this example both filter types, Baseline and RC, are compared to no filtering. The mobile phase was 1ml/min water. A length of wide ID tubing was inserted between the injection valve and the detector to simulate broad peaks. The Spray Chamber and Drift Tube temperatures were selected to provide significant noise for the example, SC 30, DT 56.


## Filter Weight



In these example the affect of the filter weight for both types is shown. Again, the Spray Chamber and Drift Tube temperatures, SC 30, DT 53, were deliberately selected to provide significant noise for the example. The mobile phase was 1ml/min water. The analyte was 100ng NaBenzoate.

This example shows the Baseline filter at filter 0 and 5

Affect of filter weight on the RC or FLT Filter type.



This example shows that while a larger filter filter weight with RC filter type decrease the noise it also broadens and decreases the peak heights.

Spray Chamber and Drift Tube temperatures, SC 30, DT 53, mobile phase was 1ml/min water analyte 1000ng NaBenzoate.