

Bende Zsolt Levente

cégvezető

06-30/311-8358 ⊠ bende.zsolt@labex.hu C 06-1/206-2455 ⊠ labex@labex.hu ⊕ www.labex.hu

LAB-EX Laborkereskedelmi Kft.

1119 Budapest, Hadak útja 41.

NextGen update

Advanced CombiFlash Techniques NextGen New Technologies

Agenda

- Sample Introduction
- Detector Options
- PeakTrak Basics
- RFID Advantages
- Chromatography Options



CombiFlash NextGen Family

	NextGen	NextGen300+	NextGen300
Flow Rate	1-100 ml/min	1-300 ml/min	1-300 ml/min
Pressure psi (bar)	150 (10.3)	300 (20)	150 (10.3) 300 (20) optional
Pumping	Single HPLC	Dual syringe	Dual syringe
Detection	200-400 nm 200-800 nm ELSD MS	200-400 nm 200-800 nm ELSD MS	200-400 nm 200-800 nm ELSD MS
Injection Valve	No	Yes	Optional
Gradients	2 solvent binary	4 solvent binary with 3 rd solvent modifier	4 solvent binary with 3 rd solvent modifier
Level Sensing	No	Yes	Yes
Air Purge	No	Yes	Yes
RFID/Level Sensing	No	Yes	Yes
User Interface	12"	12" or 15"	12"

Specs cont.

	NextGen	NextGen300+	NextGen300			
Sample Injection	Solid or liquid manual injection	Automated, self cleaning valve with injection through cartridge, syringe or directly on column	Solid or liquid manual injection. Optional automated valve			
Solvent management	None	Active level sensing on inlet and waste	Active level sensing on inlet and waste			
Dimensions (WxDxH)	14.1x17x26 in (36x43x66 cm)					
Weight	61 lb (27.7 kg) 74 lb (33.6 kg) w/ ELSD					
RFID	None	Column and racks	Optional column and racks			
Certification	CE, NRTL	CE, NRTL	CE, NRTL			
TELEDYNE ISCO Everywhereyoulook						

Sample Introduction

Good Chromatography Starts at the Beginning



Injection Parameters

	Ainimum Run Requirements Solvent Requirement: 1-Water: 0.3 L	ressary solvent requirements			
	2-Methanol: 0.4 L Expected waste: 0.7 L The run will take approximately	v 18 minutes.	Sample Loading		Solid (Pause) goes through solid load cartridge but waits – this saves time for sample prep
	Start Rack V	Left rack: 18 mm x 150 m	Solid Liquid		Solid goes through solid load cartridge and assumes sample is ready
am	Start Tube 3	Right rack: Not Detected Right r	None (on Column)		Liquid injection through top of injection valve Sample is on the
	ОК	Peak Collectio	Cancel	Ľ.,	previous run or manually loaded
olve	ont B V			37	
	TELEDYNE ISO Everywhereyou	CO Ilook			

Solid Loading

(* RealSap A

departs .

- User packed solid load cartridges
- Pre-packed solid load cartridges
- AN29 -
 - Loading Techniques
- AN88 -
 - C18 Loading Techniques
 - Celite packed cartridges
- AN82 -
 - Solid Load vs. Reservoirs





User Packed Solid Load Cartridges

- Weigh sample
- Determine amount of sorbent
 - 20% load for silica
 - 10% for Celite
- Dissolve sample completely
- Mix with silica
- Dry on rotary evaporator
- Place in cartridge
- Add top frit
- Insert solid load cartridge cap



Frit for 5 g solid load cartridge fits many 24/40 adapters commonly found in the United States and keeps silica and sample in the evaporation flask



User Packed Solid Load Cartridges

- Best resolution
- Match sample size to silica for loading
- Allows use of strong solvents
- Act as pre-column
- Prevents clogs from precipitated sample
- Also called "empty" solid load cartridges





Comparison of Sample Loading Techniques for Silica Columns



- All runs used standard solvent method on a CombiFlash
- All runs used 2% sample load
- Solvents were hexane/ethyl acetate



Pre-packed Solid Load Cartridge



Liquid Loading



Liquid loading shows reduced resolution when sample is dissolved in a strong solvent



Liquid Injection Onto Column



Sample Introduction Techniques

	Solid (Pause)	Solid	Liquid	None (Direct to Column)	
Automated		1			
Low Solubility	1	1			
Repeatable	1	1			
No carry over between runs	1	1		1	
100% sample transfer/ high recovery	1	1			
Minimal sample prep			1	1	
Tight loading band *technique dependent	* 🗸	* 🗸	* 🗸	* 🗸	2
Filters and removes baseline material	1	1			
Improved separation w longer bed	1	1			
Increases life of reusable columns (C18)	1	1			



Flash Column Load

Column Size	Easy separation Gold RediSep (Δ CV=6) :	Easy separation (∆CV=6) :	Difficult separation $(\Delta CV=1)$:
	20% loading	10% loading	1% loading
4 g	800 mg	400 mg	40 mg
12 g	2.4 g	1.2 g	120 mg
24g	4.8 g	2.4 g	240 mg
40 g	8.0 g	4.0 g	400 mg
80 g	16.0 g	8.0 g	800 mg
120 g	24 g	12 g	1.2 g
220g	44 g	22 g	2.2 g
330 g	66 g	33 g	3.3 g



Detector Options

What you see is what you get



Photo Diode Array Detector

- Variable Wavelength 200-400 nm
- Change Wavelength during the run
 - 2 wave length triggered collection
 - All wavelength triggered collection
- Purity Indicators
 - Dual Wavelength Monitoring
 - Purity Ratio



Greater Dynamic Range

- NextGen has improved UV and UV-Vis detectors
 - Up to 4 AU dynamic range
- 0.1mm fixed path length flow cell
 - Fixed creates reliable and reproducible results
 - 0.1 mm keeps compounds on scale
 - 0.3, 0.5, 1 mm optional for increased sensitivity



What You See – Detection Options

 Options allow for greatest sample recovery



All Wavelength Detection

An average of signals from a given wavelength range

Increases UV signal from low UV absorber Decreases UV signal from high UV signals to keep on scale

Average of signals



All Wavelength Collection --Application Note 81

 All-wavelength collection is a technique that monitors all detector wavelengths in a user-defined range



Purification of Green Tea Extracts Using All-Wavelength Collection

All-Wavelength Collection of "Hidden" Peaks

Purification on C18



All-Wavelength Collection Saturated Peaks

• Trace at 207 nm was saturated.

Hexane:Acetone gradient

- All-Wavelength Collection triggers collection accurately to minimize cross contamination
- Solvent Artifact is seen until the algorithm eliminates the signal





Baseline correction

- Solvent is subtracted from the base line
- See more of the compound
- Faster response time and increased recovery of compounds



Comparison Baseline Correction

Hexane/Ethyl Acetate gradient. Yellow All Wavelength detection 200-300 nm



Spectral Data





Detection Options

- Photo diode array
- Variable wavelength detection (200-400 nm) standard
- Visible wavelength detection (200-800 nm)
- All-Wavelength Detection
- ELSD integrated in unit
- Mass directed collection
 - "S" up to 1200 Daltons
 - "L" up to 2000 Daltons

Integrated ELSD Detection of non-UV Absorbing Compounds --Application Note 80

 Reaction mixture where product has no UV chromophore



ELSD Specs

	Internal ELSD
Gas Inlet Pressure	60-70 psi
Gas Consumption	<2.5 SLPM
Spray Chamber Temperature	Setting range 10-60 °C
Drift Tube Temperate	Setting range 30-90 °C
Split Flow Rate	0.75 ml/min 5-30 ml/min Increases from 0.75 to 1.5 ml/min 30-60 ml/min 1.5 ml/min 60-300 ml/min
Sensitivity	High & Normal settings

ELSD detection is not recommended for flow rates below 5 ml/min

ELSD can be turned off for any runs below this recommendation. -Note flash columns \geq 4 grams are optimal over 15 ml/min.



New ELSD Spray Chamber



- New spray chamber compatible with aqueous buffers
- Design is simple to clean
- Offers <u>High</u> and <u>Normal</u> sensitivity software selection

Integrated ELSD available on all CombiFlash models

- No additional lab space required
- Available on NextGen 300+, 300, 100



Two Modes of Operation

- Normal Sensitivity
 - Set to 1/2 of Lumen sensitivity
 - Gain default set to 2X
 - Results in same peak magnitude as current response
 - Ability to lower the sensitivity for strong signals
- High Sensitivity
 - Set to 6X Lumen sensitivity
 - Ability to cut peaks of <5mg samples



ELSD Sensitivity Selection

Normal Sensitivity ELSD Signal

High Sensitivity ELSD Signal



Improved Recovery over UV Detection

 Shaded areas on peaks denote peak area detected by ELSD but not with UV



Improved Recovery over Competitors

A sample of 2,3isopropylidene-Dribofuranose was synthesized and the reaction mixture (1.0 mL) was purified on a 40g RediSep Rf silica column (PN 69-2203-340) using the default gradient on each flash system. Gradient was hexanes/ethyl acetate.





CombiFlash PurIon MS

- Available on all Combiflash NextGen, EZPrep, ACCQPrep models
- A fully integrated operation in PeakTrak software
 - Simple to operate
- Mass directed purification
- Compound verification
- UV/Vis and ELSD detection along with MS





Purlon Module Specs

	"S" Model	"L" Model		
Upgrade available	Rf+, Lumen, EZPrep, ACCQPrep, NextGen			
Dimensions (WxHxD)	11x26x22 in (28x66x56 cm) MS and FI 9x10x18 in (23x26x46 cm) Roughing Pump			
MS Detection	10-1200 Dalton 1 Dalton resolution	10-2000 Dalton 1 Dalton resolution		
Ionization Mode	Positive and Negative Automatic Switching			
Capillary Cone Removal	Does not require vacuum shutdown			
Probes	ESI or APCI			
Sample injection	Direct Injection for MS Method	Development		



Benefits to the End User

- Determine presence of co-eluting impurity
 - Quickly determine that additional purification needed
- Collect only on target mass
 - Increase throughput
 - Shorten run times by terminating purification after target has eluted

37

- Lower solvent consumption and costs
- Eliminate the need for post-purification analysis
 - No more TLC
 - No post-run LC/MS
- Collect over a mass range
 - Eliminate unwanted masses
 - Useful in Natural Product Research



Terminate On Target

- Use MS Method Development screen
- Select desired masses





Sample run

- Weak UV absorbance
- Run entire program



Terminate on Target

- Terminate run after peak with m/z=180 Da collected
- Collection with UV and MS
- Time and solvent saved



Terminate on Target

- Second peak collected, terminate on m/z = 185
- Collection on SIM (single ion monitoring) 185 Da



Terminate on target

- Terminate run after m/z 180 AND 185 Da collected.
- Up to 4 SIM or 3 SIM and 1 XIC (extracted ion current)



What You Get – Optimized Fraction Collection







Peak Parameters

- Slope Based
 - Peaks will be recorded if the slope algorithm indicates a peak within the Peak Width setting
- Peak Width
 - Average peaks are detected based on a range from 0.2 to 2 times the peak width at baseline.
- Threshold
 - Collects peaks if the absorbance unit (AU) set is exceeded.
- Select both Slope and Threshold as an OR operation
- Monitor
 - Detection signal will be displayed but will not trigger collection



PeakTrak Basics



Basic PeakTrak Run Screen



PeakTrak Main Screen



PeakTrak Method Editor

Flow rate – established in the column method





PeakTrak Method Editor



Behind the scenes-Configuration of CombiFlash

- Select solvent per inlet line
- Add and remove solvents
- Enable solvent sensing (300+)
- Set solvent detect level (300+)
- Turn on audible alarms
- Name the CombiFlash
- Set the time zone
- Set vapor limit
 - high/med/low/off
 - optional
- Set default tube volume





Network Configuration

- View the unit from desktop
- Print via network printer
- Auto save to network
- Need to work with local IT department
- TN 28 Networking Guidelines

Configuration	
Instrument IP address configuration	Network V DHCP
Enter the URL for the printer below. A test page will be Leave the URL blank to turn off network printing. For a LPR print queue use the form: For a LPR print queue use the form: Ipd://ip_address/queue_name socket Printer Type Connection URL PostScript Connection URL	e printed when you save these settings. JetDirect print queue use the form: et://ip_address:port
Network file save configuration File Type Domain Disabled tdv.teledvne.com Network Share //lin-vpsvs-wfs06/rspatz\$/netSave OK	Admin Server User Name Spatz Cancel



User Preferences

- Set Language
- Auto print function
- Default display units
- Enable peak hold to extend the run
- Enable rapid equilibration
- Enable auto peak holds
- Pressure in psi or bar

	Configuration
	Instrument Configuration Network Configuration User Preferences
	Changing preferences for user: Common
	LanguageAutomatically printEnglishreport at end of run
d	Default Run Units Time, Minutes Enable run length extension
	Pressure Units PSI DEnable rapid equilibration
	Enable automatic peak hold
	OK Cancel



Post Run – View Separations



View Separations



Help Window

Method Editor	Tools	Help		
od Name: Temporary 5		Service		
		Software Update		
		Expor	t Log Files	
			eorder Code	
		Abou	: PeakTrak	
		Abou	: Teledyne ISCO	



Everywhereyoulook"



- Service access the service functions
- Software Update –update from a USB drive
- Export Log Files export to troubleshoot any software bugs
- Last Reorder Code active for Column Plan (300+)
- About PeakTrak tells the software version
- About Isco Contact Teledyne Isco

RFID Technology



RediSep Information



Scale Up from Milligrams to Grams

- Easily convert 4 g column method to 330 g column method on the Rf and up to 3 kg on Torrent
- Open method to be scaled
- Go to "File" then "Save method as..."

FILE METHOD EDIT	TOR TOOLS HELP	FILES				
0.50		/ comr	non /	COPY PASTE DELETE		
0.45			E NAME		E	
0.40			SubFolder1			
0.30		M	base.mtd	2010-04-2	3	
0.25			test1 mtd	2010-02-1	2	
0.20			Burnla actal	2010-02-1	2	
0.10			Burnin.mta	2010-02-1	2	
0.05						
0.00			Silica 80g	Silica		• • •
0.0 ABSORBANCE	5.0		Silica 120g	g Silica		25.0 TH
C	ZOOM IN GRAPH	POINT	Silica 220	j Silica		Min.
REDISEP			Silica 330ç	g Silica		
Silica 120g			Silica 750ç	g Silica		
FLOW RATE 85 ml/mir	Wavelength	1 (254nm)	Silica 1.5k	g Silica	SAVE	
SOLVENT A	SOLVENT B		SCALE NEW REDISE	P SIZE	CANCEL	NONE
hexane	ethyl acetate		Silica 120			

4g column



Scaled to 330g Column



Scaled to 1.5 kg Column



Load Method from Column



Load Method from Rack



Chromatography Options



LC Comparisons

LC Method	Teledyne ISCO	Particle Size	Flow Rate Range	Back Pressure	Capacity
Process Flash	Torrent/ NextGen	20-60 μ	0.3-1 L/min	25-75 PSI	g-100's g
Flash	NextGen/ EZPrep	20-60 μ	>18 mL/min	10-50 PSI	mg-g
Prep HPLC	EZPrep/ ACCQPrep	5-20 μ	3-100 mL/min	200 – 2000 PSI	10-500 mg
Analytical HPLC	NA	5 u	1 mL/min	500- 6000 PSI	<5 mg
UHPLC	NA	1.3-3 μ	0.1-0.5 mL/min	>6000 PSI	



Thank you!!

Guidelines & Tactics for Flash Chromatography

Free download www.teledyneisco.com/flashguide

