



## -86 °C ULT FREEZER

Series	MODEL	Type	Insulation Type	High Energy	Cooling system	Consumables	REMARK
						Rack Quantity	
Energy saving series	DW-86L486E	Upright	VIP	Yes	Cascade	DCJ-54-A(B) * 16	HC gas, keyboard controller
	DW-86L628E	Upright	VIP	Yes	Cascade	DCJ-54-A(B) * 20	HC gas, keyboard controller
	DW-86L338J	Upright	VIP	Yes	Cascade	DCJ-54-A(B) * 6 + DCJ-44-A(B) *6	HC gas, keyboard controller
	DW-86L388J	Upright	VIP	Yes	Cascade	DCJ-55-A(B) * 12	HC gas, keyboard controller
	DW-86L490J	Upright	VIP	Yes	Cascade	DCJ-54-A(B) * 8 + DCJ-44-A(B) *8	HC gas, keyboard controller, double door
	DW-86L578J	Upright	VIP	EEC	Cascade	DCJ-55-A(B) * 16	HC gas, keyboard controller
	DW-86L728J	Upright	VIP	EEC	Cascade	DCJ-55-A(B) * 20	HC gas, keyboard controller
	DW-86L828J	Upright	VIP	EEC	Cascade	DCJ-55-A(B) * 24	HC gas, keyboard controller
	DW-86L100J	Under Bench	VIP	Yes	Cascade	DCJ-53-A(B) * 4	HC gas, keyboard controller
	DW-86W100J	Chest	Standard	Yes	Cascade	DCJ-08-A * 9	HC gas, keyboard controller
DW-86W420J	Chest	Standard	Yes	Cascade	DCJ-10-A * 27	HC gas, keyboard controller	
Water cooling Series (Touch LCE)	DW-86L828W	Upright	VIP	Yes	Water cooling	DCJ-55-A(B) * 24	HC gas, touch screen controller, inverter technology
Water cooling Series (Touch LCE)	DW-86L959W	Upright	VIP	Yes	Water cooling	DCJ-55-A(B) * 28	HC gas, touch screen controller, inverter technology
Double Cooling Series	DW-86L578S	Upright	VIP	EEC	Double Cooling	DCJ-55-A(B) * 16	HC gas, keyboard controller
	DW-86L728S	Upright	VIP	EEC	Double Cooling	DCJ-55-A(B) * 20	HC gas, keyboard controller
	DW-86L578ST	Upright	VIP	EEC	Double Cooling	DCJ-55-A(B) * 16	HC gas, touch screen controller
	DW-86L728ST	Upright	VIP	EEC	Double Cooling	DCJ-55-A(B) * 20	HC gas, touch screen controller
Smart Frequency Series (LED)	DW-86L579BP	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 16	HC gas, keyboard controller

<b>Smart Frequency Series (LED)</b>	DW-86L729BP	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 20	HC gas, keyboard controller
<b>Smart Frequency Series (LED)</b>	DW-86L829BP	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 24	HC gas, keyboard controller
	DW-86L959BP	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 28	HC gas, keyboard controller
<b>Smart Frequency Series (Touch LCD)</b>	DW-86L579BPT	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 16	HC gas, touch screen controller
	DW-86L729BPT	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 20	HC gas, touch screen controller
	DW-86L829BPT	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 24	HC gas, touch screen controller
	DW-86L959BPT	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 28	HC gas, touch screen controller



**LAB-EX Laborkereskedelmi Kft.**

1119 Budapest, Hadak útja 41.

**Tel.:** 06-1-206-2455

**E-mail:** labex@labex.hu

**Web:** www.labex.hu

