












## Rendkívül kedvezményes árak az alábbi DEMO pumpákra:

SSI stock pumps	Description	Price (without VAT)	Product details	
LD012PFP3A	LD Class Dual Piston Pump with Belt Drive and Shaped Cams; 0.000 - 12.000 mL/min.; 6000 psi; PEEK Fluid Path; Automatic Piston Wash; with Pulse Dampener; Standard Seal (Flouropolymer energizer); Inlet (1/8" OD Flexible Tubing); Outlet	1.650.000,-Ft	<a href="https://www.teledyneisco.com/en-us/pumps/reaxus-ld-class">https://www.teledyneisco.com/en-us/pumps/reaxus-ld-class</a>	
LS010SFT2ACT	LS Class Single Piston Pump with Belt Drive and Shaped Cam; 0.00 - 10.00 mL/min.; 6000 psi; Stainless Steel Fluid Path; Automatic Piston Wash; with Pulse Dampener; Organic Seal (stainless steel energizer); Inlet (1/8" OD Flexible Tubing); Outlet (1/16" OD Tubing)	1.163.000,-Ft	<a href="https://www.teledyneisco.com/en-us/pumps/reaxus-ls-class">https://www.teledyneisco.com/en-us/pumps/reaxus-ls-class</a>	
M1010SFN2ACT	ChromTech M1 Class 10ML SS SF NPD ORS FB-10 OEM pump	1.140.000,-Ft	<a href="https://www.teledyneisco.com/en-us/pumps/reaxus-m1-class">https://www.teledyneisco.com/en-us/pumps/reaxus-m1-class</a>	
MX010SFP4C	MX Class Single Piston Pump with Eccentric Drive and Shaped Cam; 0.00 - 10.00 mL/min.; 5000 psi; Stainless Steel Fluid Path; with automatic Piston Wash; with Pulse Dampener Only; Organic Seal (flouropolymer energizer); Inlet (1/8" OD Rigid Tubing customer supplied); Outlet (1/8" OD Tubing customer supplied)	1.000.000,-Ft	<a href="https://ssihplc.com/product-category/hplc-pumps/mx-class/">https://ssihplc.com/product-category/hplc-pumps/mx-class/</a>	

MX010SFT1A	MX Class Single Piston Pump with Eccentric Drive and Shaped Cam; 0.00 - 10.00 mL/min.; 5000 psi; Stainless Steel Fluid Path; Automatic Piston Wash; with Pulse Dampener; Standard Seal (stainless steel energizer); Inlet (1/8" OD Flexible Tubing); Outlet (1/16" OD Tubing)	1.000.000,-Ft	<a href="https://ssihplc.com/product/mx010sft1a/">https://ssihplc.com/product/mx010sft1a/</a>	
LD024SFX1A	LD Class Dual Piston Pump with Belt Drive and Shaped Cams; 0.000 - 24.000 mL/min.; 10,000 psi; Stainless Steel Fluid Path; Automatic Piston Wash; with Pressure & Leak Sensors, no dampener; Standard Seal (stainless steel energizer); Inlet (1/8" OD Flexible Tubing); Outlet (1/16" OD Tubing)	1.650.000,-Ft	<a href="https://www.teledyneisco.com/en-us/pumps/reaxus-ld-class">https://www.teledyneisco.com/en-us/pumps/reaxus-ld-class</a>	
LD024SFT1A	LD Class Dual Piston Pump with Belt Drive and Shaped Cams; 0.00 - 24.00 mL/min.; 10000 psi; Stainless Steel Fluid Path; with Piston Wash; with Pressure & Leak Sensors, with pulse dampener; Standard Seal (stainless steel energizer); Inlet (1/8" OD Flexible Tubing); Outlet (1/16" OD Tubing)	1.680.000,-Ft	<a href="https://www.teledyneisco.com/en-us/pumps/reaxus-ld-class">https://www.teledyneisco.com/en-us/pumps/reaxus-ld-class</a>	
S10SNXP1ST	Supercritical 24 Pump, constant pressure, self-flush, 0.01-24 mL/min., 110/220V, with RS232 control, pressure monitoring, thermo-electric head chillers	2.350.000,-Ft	<a href="https://ssihplc.com/product/s10snxp1/">https://ssihplc.com/product/s10snxp1/</a>	
P25SFX01	HF Series Pump, 300mL, Stainless Steel, Self-Flush, with Pulse Dampener and Transducer, Constant Flow, Large Motor	3.600.000,-Ft	<a href="https://ssihplc.com/product/p25sfx01/">https://ssihplc.com/product/p25sfx01/</a>	

PR100SFX1A	PR Class Dual Piston Pump with Belt Drive and Shaped Cams; 0.0 - 100.0 mL/min.; 4000 psi; Stainless Steel Fluid Path; Automatic Piston Wash; with Pressure & Leak Sensors, no dampener; Standard Seal (stainless steel energizer); Inlet (1/8" OD Flexible Tubing); Outlet (1/16" OD Tubing)	1.900.000,-Ft	<a href="https://ssihplc.com/product/pr100sfx1a/">https://ssihplc.com/product/pr100sfx1a/</a>	
LS010SFT1XPTLSC E	LS Class Single Piston Pump with Belt Drive and Shaped Cam; 0.00 - 10.00 mL/min.; 6000 psi; Stainless Steel Fluid Path; Automatic Piston Wash; with Pulse Dampener; Standard Seal (stainless steel energizer); Inlet (1/8" OD Flexible Tubing); Outlet (1/16" OD Tubing)	1.163.000,-Ft	<a href="https://www.teledyneisco.com/en-us/pumps/reaxus-ls-class">https://www.teledyneisco.com/en-us/pumps/reaxus-ls-class</a>	

***Az ajánlat csak a készlet erejéig érvényes! A feltüntetett árak az áfát nem tartalmazzák.***

Kérdés, árajánlatkérés esetén vegye fel a kapcsolatot velünk!

**Kapcsolattartó:** Bende Zsolt  
**E-mail:** [bende.zsolt@labex.hu](mailto:bende.zsolt@labex.hu)  
**Mobil:** 06-30-311-8358

