

Product number	Product type	Serial Number
622130009R	FOXY R2 FC SA 100-240V REFR	215E20051
622130101R	FOXY R2 SA FC OEM 100-240 VOLTS REFR	218J03052
622130006R	FOXY R2 FC 100-240 VOLTS REFR	218J03052
623870012R	FOXY R1 SA FC/UNIT 100-240V REFR	219D02271
623870071R	FOXY R1 SA FC OEM 100-240 VOLTS REFR	219K01651

	Fraction Collectors	PRICE in USD (without VAT)
68-2130-009	Foxy R2 Fraction Collector Operating temperature range from 0 to 40 degrees C. RFID sensors for fast and error free setup. Standard 1/8 inch OD (0.062 inch ID) tubing included for collection up to 125 mL/min. Includes (2) racks for (288) 13 mm tubes. (Substitution of same-price racks must be specified on purchase order. Tubing package for 0 to 25 mL/min collection, extra.)	\$5.000
68-3870-012	Foxy R1 Fraction Collector Operating temperature from 0 to 40 degrees C. Standard 1/16 inch OD (0.030 inch ID) tubing included for collection up to 25 mL/min. Includes (1) rack for (144) 13 mm tubes. (Substitution of same-price racks must be specified on purchase order. Tubing package for 25-125 mL/min collection, extra.)	\$3.500

	Foxy® R1	Foxy® R2
Dimensions (H x W x D)	14 x 12.25 x 13 in (35.5 x 31.1 x 33.0 cm)	14.9 x 12.25 x 21 in (37.8 x 31.1 x 53.3 cm)
Weight	15.6 lbs (7.1 kg) -- with drip tray, no rack	22.8 lbs (10.3 kg) -- with drip tray, no rack
Maximum Flow Rate	25 mL/min standard; 125 mL/min with optional accessory package	125 mL/min standard; 0 - 25 mL/min with optional accessory package; 1000 mL/min with optional high flow diverter valve
Operating Voltage	100 - 240 VAC, 50-60 Hz, 1 ampere maximum	100 - 240 VAC, 50-60 Hz, 1 ampere maximum
	Brochure	

Kapcsolattartó: Bende Zsolt
Mobil: 06-30-311-8358

E-mail: bende.zsolt@labex.hu

LAB-EX Laborkereskedelmi Kft. Tel.: 06-1-206-2455 E-mail: labex@labex.hu Web: www.labex.hu