



## -86 °C ULT FREEZER

Series	MODEL	Capacity (L)	Type	Insulation Type	High Energy	Cooling system	Consumables	REMARK
							Rack Quantity	
Energy saving series	<a href="#">DW-86L338J</a>	338	Upright	VIP	Yes	Cascade	DCJ-54-A(B) * 6 + DCJ-44-A(B) *6	HC gas, keyboard controller
	<a href="#">DW-86L388J</a>	388	Upright	VIP	Yes	Cascade	DCJ-55-A(B) * 12	HC gas, keyboard controller
	<a href="#">DW-86L490J</a>	490	Upright	VIP	Yes	Cascade	DCJ-54-A(B) * 8 + DCJ-44-A(B) *8	HC gas, keyboard controller, double door
	<a href="#">DW-86L578J</a>	578	Upright	VIP	EEC	Cascade	DCJ-55-A(B) * 16	HC gas, keyboard controller
	<a href="#">DW-86L728J</a>	728	Upright	VIP	EEC	Cascade	DCJ-55-A(B) * 20	HC gas, keyboard controller
	<a href="#">DW-86L828J</a>	828	Upright	VIP	EEC	Cascade	DCJ-55-A(B) * 24	HC gas, keyboard controller
	<a href="#">DW-86L100J</a>	100	Under Bench	VIP	Yes	Cascade	DCJ-53-A(B) * 4	HC gas, keyboard controller
	<a href="#">DW-86W100J</a>	100	Chest	Standard	Yes	Cascade	DCJ-08-A * 9	HC gas, keyboard controller
<a href="#">DW-86W420J</a>	120	Chest	Standard	Yes	Cascade	DCJ-10-A * 27	HC gas, keyboard controller	
Water cooling Series (Touch LCE)	DW-86L828W	828	Upright	VIP	Yes	Water cooling	DCJ-55-A(B) * 24	HC gas, touch screen controller, inverter technology
Water cooling Series (Touch LCE)	DW-86L959W	959	Upright	VIP	Yes	Water cooling	DCJ-55-A(B) * 28	HC gas, touch screen controller, inverter technology
Double Cooling Series	<a href="#">DW-86L578S</a>	578	Upright	VIP	EEC	Double Cooling	DCJ-55-A(B) * 16	HC gas, keyboard controller
	<a href="#">DW-86L728S</a>	728	Upright	VIP	EEC	Double Cooling	DCJ-55-A(B) * 20	HC gas, keyboard controller
	<a href="#">DW-86L578ST</a>	578	Upright	VIP	EEC	Double Cooling	DCJ-55-A(B) * 16	HC gas, touch screen controller
	<a href="#">DW-86L728ST</a>	728	Upright	VIP	EEC	Double Cooling	DCJ-55-A(B) * 20	HC gas, touch screen controller
Smart Frequency Series (LED)	<a href="#">DW-86L579BP</a>	579	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 16	HC gas, keyboard controller

<b>Smart Frequency Series (LED)</b>	<a href="#">DW-86L729BP</a>	729	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 20	HC gas, keyboard controller
<b>Smart Frequency Series (LED)</b>	<a href="#">DW-86L829BP</a>	829	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 24	HC gas, keyboard controller
	<a href="#">DW-86L959BP</a>	959	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 28	HC gas, keyboard controller
<b>Smart Frequency Series (Touch LCD)</b>	<a href="#">DW-86L579BPT</a>	579	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 16	HC gas, touch screen controller
	<a href="#">DW-86L729BPT</a>	729	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 20	HC gas, touch screen controller
	<a href="#">DW-86L829BPT</a>	829	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 24	HC gas, touch screen controller
	DW-86L959BPT	959	Upright	VIP	EEC	variable speed drive control	DCJ-55-A(B) * 28	HC gas, touch screen controller



**LAB-EX Laborkereskedelmi Kft.**  
1119 Budapest, Hadak útja 41.

**Tel.:** 06-1-206-2455    **E-mail:** labex@labex.hu    **Web:** www.labex.hu

