

Which is the most suitable RediSep® Column for me?

	Gold	Silver	Bronze
Media	Silica	Silica	Silica
# Runs	Single use	Single use	Single use
Size range	4 g–3 kg	4 g–3 kg	4–330 g
Particle Size	20–40 µm	40–63 µm	40–60 µm
Mesh size	400–632	230–400	230–400
Pore Size	60 Å	60 Å	55–60 Å
pH	6–8	6–8	6.7
Particle shape	Spherical	Irregular	Irregular
Loading capacity	0.1–20%	0.1–10%	0.1–8%
Pressure Rating* up to 300 g	600–210 psi	600–210 psi	600–210 psi
RFID (Teledyne ISCO systems)	Yes	Yes	No
Manufacture Location	USA	USA	USA
QA testing	5%	1%	0.1%
Warranty—any damaged column we ship a free box of columns	Yes	Yes	No
Quantity	X	XX	XXX
Price	\$\$\$	\$\$	\$

Kapcsolattartó: Bende Zsolt

E-mail: bende.zsolt@labex.hu

Mobil: 06-30-311-8358