

# XYSCALE

WWW.XYSCALE.COM



WaterProof Scale



# ÁRLISTA 2021







Analytical Balance











Moisture Analyzer









Weighing Scale

NAME	PICTURE	MODE	PARAMETERS	Price in HuF (without VAT)	Pan size	Note
<b>SEMIMICRO BALANCES</b>						
<b>FA-SEM</b> Series Semi micro Balance		FA1265SEM	62g/120g 0.01mg/0.1mg	509 286 Ft	Φ80mm	1. Electromagnetic balanced sensors 2. Automatic internal calibration 3. Dimension: 345×215×345mm 4. Efficient filtering to deliver quick and accurate results 5. LCD with backlight with adjustable contrast 6. Glass draft shield with 3 sliding doors for easy access to the items being weighed 7. Plastic protection cover for keypad 8. -ION: with lonser
		FA2285SEM	82g/220g 0.01mg/0.1mg	553 571 Ft		
		FA22105SEM	102g/220g 0.01mg/0.1mg	620 000 Ft		
		FA125SEM-ION	122g/0.01mg	620 000 Ft		
		FA165SEM-ION	162g/0.01mg	664 286 Ft		
		FA205SEM-ION	202g/0.01mg	752 857 Ft		
<b>ANALYTICAL BALANCE</b>						
<b>FA-AI</b> Series Analytical Balance		FA214AI	220g/0.1mg	332 143 Ft	Φ80mm	1. Electromagnetic balanced sensors 2. Automatic internal calibration 3. Dimension: 345×215×345mm 4. Efficient filtering to deliver quick and accurate results 5. LCD with backlight with adjustable contrast 6. Glass draft shield with 3 sliding doors for easy access to the items being weighed 7. Plastic protection cover for keypad
		FA314AI	310g/0.1mg	343 214 Ft		
		FA414AI	410g/0.1mg	442 857 Ft		
		FA514AI	510g/0.1mg	487 143 Ft		
		FA614AI	610g/0.1mg	531 429 Ft		
<b>ANALYTICAL BALANCE</b>						
<b>FA-T</b> Series Analytical Balance		FA1004T	100g/0.1mg	148 357 Ft	Φ80mm	1. Electromagnetic balanced sensors; 2. Tare, counting, unit conversion, under weighing functions; 3. RS232 interface (optional). 4. <b>5 inch</b> touch panel 5. Density test/Counting Test/Percentage Test 6. <b>Internal Calibration</b>
		FA2004T	200g/0.1mg	159 429 Ft		
		FA2204T	220g/0.1mg	168 286 Ft		
<b>FA</b> Series Analytical scales		FA1004	100g/0.1mg	104 071 Ft	Φ80mm	1. Electromagnetic balanced sensors; 2. Tare, counting, unit conversion, under weighing functions; 3. RS232 interface (optional). 4. External Calibration
		FA1204	120g/0.1mg	106 286 Ft		
		FA2004	200g/0.1mg	110 714 Ft		
		FA2204	220g/0.1mg	115 143 Ft		
		FA1004N	100g/0.1mg	126 214 Ft	Φ80mm	1. Electromagnetic balanced sensors; 2. Tare, counting, unit conversion, under weighing functions; 3. RS232 interface (optional). 4. <b>Internal Calibration</b>
		FA1204N	120g/0.1mg	126 214 Ft		
		FA2004N	200g/0.1mg	128 429 Ft		
		FA2204N	220g/0.1mg	132 857 Ft		
<b>ANALYTICAL BALANCE</b>						
<b>JA</b>		JA203	200g/1mg	101 857 Ft	Φ80mm	1. Electromagnetic balanced

Series Analytical scales		JA303	300g/1mg	97 429 Ft		sensors; 2. Tare, counting, unit conversion, under weighing functions; 3. RS232 interface (optional). 4.External Calibration
		JA503	500g/1mg	104 071 Ft		
		JA603	600g/1mg	106 286 Ft		
		<b>JA1003</b>	<b>1000g/1mg</b>	139 500 Ft		
		<b>JA1203</b>	<b>1200g/1mg</b>	155 000 Ft		
		JA203N	200g/1mg	119 571 Ft		
		JA303N	300g/1mg	119 571 Ft		
		JA503N	500g/1mg	126 214 Ft		
		JA603N	600g/1mg	130 643 Ft		
		JA103H	110g/1mg	55 357 Ft		Φ90mm
JA203H		210g/1mg	57 571 Ft			
JA303H		310g/1mg	62 000 Ft			
<b>ANALYTICAL BALANCE (Economy Mode)</b>						
FA-E JA-E Economy Balance		FA1004E	100g/0.1mg	73 071 Ft	Φ80mm Option: Φ90mm	*5 glass draft shield *Count/Percentage/Check/Peak *RS232 interface(Optional USB) *Option internal cal(N)+300 *Option under weighing+100 *Option density function (MOQ:10PCS)
		FA1104E	110g/0.1mg	73 071 Ft		
		FA1204E	120g/0.1mg	73 071 Ft		
		FA2004E	200g/0.1mg	77 500 Ft		
		FA2104E	210g/0.1mg	78 829 Ft		
		FA2204E	220g/0.1mg	81 043 Ft		
		JA103E	100g/1mg	73 071 Ft		
		JA203E	200g/1mg	73 071 Ft		
		JA303E	300g/1mg	77 500 Ft		
		JA503E	500g/1mg	77 500 Ft		
<b>PRECISION BALANCE</b>						
JA-P Electronic Balance		JA103P	110g/0.001g	32 329 Ft	Φ90mm	1.Unit conversion; 2. LCD digital display; 3.Tare, counting, under weighing functions; 4.RS232 interface (optional). 5.Rechargeable battery(optional) 6.Option dry battery and lithium battery
		JA203P	210g/0.001g	34 543 Ft		
		JA303P	310g/0.001g	37 643 Ft		
		JA103PL	100g/0.001g 600g/0.01g	34 543 Ft		
<b>PRECISION BALANCE</b>						
XY-C Series Precision scales		XY200C	210g/0.01g	13 729 Ft	Φ130mm	1.Stainless steel Pan, metal coating base; 2. LED digital display; 3.Tare, counting functions; 4.RS232 interface (optional); 5. battery (optional); 6. windbreak mask (optional). 10pcs in one box
		XY300C	310g/0.01g	13 729 Ft		
		XY600C	610g/0.01g	13 729 Ft		
		XY1000C	1100g/0.01g	15 500 Ft		
		XY2000C	2100g/0.01g	23 914 Ft		
		XY3000C	3100g/0.01g	25 686 Ft		
<b>PRECISION BALANCE</b>						

<b>XY-2C</b> <b>Series</b> <b>Precision</b> <b>scales</b>		XY2002C	210g/0.01g	13 729 Ft	Φ130mm	1.Stainless steel Pan; 2.LCD digital display; 3.Tare, counting functions; 4.RS232 interface (optional); 5.battery (optional); 6. windbreak mask (optional). unit: oz,ct,lb,g 10pcs in one box
		XY3002C	310g/0.01g	13 729 Ft		
		XY6002C	610g/0.01g	13 729 Ft		
		XY10002C	1100g/0.01g	15 500 Ft		
		XY20002C	2100g/0.01g	23 914 Ft		
		XY30002C	3100g/0.01g	25 686 Ft		
<b>JA-C</b> <b>Series</b> <b>Precision</b> <b>scales</b>		JA3000C	3100g/0.01g	57 571 Ft	180*180mm	1.Aluminium shell, anti-jamming, anti-radiation; 2. LED digital display; 3.Tare, counting functions; 4.RS232 interface (optional);
		JA4000C	4100g/0.01g	62 000 Ft		
		JA5000C	5100g/0.01g	62 000 Ft		
<b>PRECISION BALANCE</b>						
<b>XY-BF</b> <b>Series</b> <b>Precision</b> <b>scales</b>		XY2000B	2100g/0.1g	15 500 Ft	Φ130mm	1.Stainless steel Pan, metal coating base; 2. LED digital display; 3.Tare, counting functions; 4. RS232 interface (optional); 5. battery (optional); 6. windbreak mask (optional). 10pcs in one box
		XY3000BF	3100g/0.1g	15 500 Ft	160*160mm	
		XY5000BF	5100g/0.1g	15 500 Ft		
		XY10KBF	11kg/0.2g	24 357 Ft		
<b>XY-1B</b> <b>Series</b> <b>Precision</b> <b>scales</b>		XY6001B	610g/0.1g	10 186 Ft	Φ130mm	1. Stainless steel Pan; 2.LCD digital display; 3.Tare, counting functions; 4. RS232 interface (optional); 5. battery (optional); 6.windbreak mask (optional). 4unit: oz,ct,lb,g 10pcs in one box
		XY10001B	1100g/0.1g	10 629 Ft		
		XY20001B	2100g/0.1g	15 500 Ft		
		XY30001BF	3100g/0.1g	15 500 Ft	160*160mm	
		XY50001BF	5100g/0.1g	15 500 Ft		
<b>MOISTURE ANALYZER</b>						
<b>XY-MW</b> <b>Halogen</b> <b>Moisture</b> <b>analyzer</b>		XY-110MW	110g/10mg	146 143 Ft	Φ90mm	1.Halogen heating 2.Pansize:90mm 3.Stainless steel chamber 4.Temperature and time can be set 5.The percentage of moisture content 6.40-199℃ 7.XY-100MW-IR:Infrared lamp
		<b>XY-105MW</b>	110g/5mg	165 200 Ft		
		XY-102MW	110g/2mg	194 857 Ft		
		<b>XY-100MW</b>	110g/1mg	240 000 Ft		
		XY-100MW-IR	110g/1mg	256 857 Ft		
<b>XY-MW-T</b> <b>Halogen</b> <b>Moisture</b> <b>analyzer</b>		XY-100MW-T	110g/1mg	310 000 Ft	Φ90mm	1.Halogen heating 2.Pansize:90mm 3.Stainless steel chamber 4.Temperature and time can be set 5.The percentage of moisture content 6. 5inch touch panel 7.Analysis curve 8..40-199℃
<b>XY-MS</b> <b>Halogen</b> <b>Moisture</b> <b>analyzer</b>		XY-100MS	110g/1mg	487 143 Ft	Φ100mm	1.Halogen heating 2.Pansize:100mm 3.Stainless steel chamber 4.Temperature and time can be set 5.The percentage of moisture content 6.7inch touch screen 7.Automatic open... 8.Easy cleaning 9.230degree heating temperature... Plus( Rechargeable battery+wireless)
		XY-100MS (Plus)	110g/1mg	522 571 Ft		
		XY-200MS	210g/1mg	553 571 Ft		
		XY-200MS (Plus)	210g/1mg	597 857 Ft		
<b>HYDROSTATIC WEIGHING SCALE</b>						

XY- J Hydrostatics balance		XY2000J	2100g/0.1g	37 643 Ft	160*160mm	1.Stainless steel Pan, metal coating base 2.LED digital display; 3.Tare, counting functions; 4.Samples density test.
		XY3000J	3100g/0.1g	37 643 Ft		
		XY5000J	5100g/0.1g	42 071 Ft		
<b>WEIGHING SCALE</b>						
XY- M Series Precision scales		XY10MA	11kg/1g	27 457 Ft	315*230mm	1.Stainless steel big Pan; 2.LCD digital display; 3.to have skin, counting functions; 4.RS232 interface (optional); 5.Rechargeable battery (optional). 4pcs in one box
		XY15MA	16kg/1g	27 457 Ft		
		XY20MA	21kg/1g	31 000 Ft		
		XY30MA	31kg/1g	31 000 Ft		
		XY5MB	5.1kg/0.1g	25 686 Ft		
		XY10MB	11kg/0.1g	31 000 Ft		
		XY15MB	16kg/0.1g	34 100 Ft		
		XY20MB	21kg/0.1g	34 100 Ft		
		XY25MB	25kg/0.1g	41 186 Ft		
		XY28MB	28kg/0.1g	44 286 Ft		
XY- EF Series Precision scales		XY60E	61kg/10g	34 100 Ft	300×400 410×510 450×600 600×800 mm	1.Stainless steel Pan, plastic spraying carbon steel base big; 2. LED digital display; 3. Tare, counting functions; 4. RS232 interface (optional); 5. Rechargeable Battery (optional). 6.410X510 standard pan size
		XY100E	110kg/10g	34 100 Ft		
		XY150E	160kg/10g	34 100 Ft		
		XY200E	210kg/10g	37 643 Ft	600×800mm	
		XY300E	310kg/10g	37 643 Ft		
		XY500E	510kg/50g	70 857 Ft		
		XY60F	61kg/1g	42 071 Ft	410×510mm	
		XY100F	110kg/1g	44 286 Ft		
		XY150F	160kg/1g	55 357 Ft		
<b>WEIGHING SCALE</b>						
XY- MT Internal Weighing Scale		XY10MTA	10kg/1g	93 000 Ft	305x230 mm	1. 7 Inch touch panle 2.*Weighing,density,check weight,label print,%,counting 3.Internal calibration(external reduce \$50 USD) 4.USB input and output 5.Information enter 6.Internal calibration 【≤10kg 0.1g/30kg 1g】 * ... ..
		XY30MTA	30kg/1g	97 429 Ft		
		XY10MTB	10kg/0.1g	97 429 Ft		
		XY20MTB	20kg/0.1g	104 071 Ft		
		XY30MTB	30kg/0.1g	132 857 Ft		
		XY3MTC	3kg/0.01g	146 143 Ft		
		XY4MTC	4kg/0.01g	155 000 Ft		
		XY5MTC	5kg/0.01g	159 429 Ft		
<b>WATERPROOF WEIGHING SCALE</b>						
XY- WP Waterproof Scale		XY3KB1WP	3kg/0.1g	46 500 Ft	305x230 mm	1. Stainless steel bracket 2. Waterproof box+Gluing protect display 3. Sealing glue protect main board 4. Dual display 5. 0.4s stable time 6. Washable IP68/CE/RoHS/LFGB certification
		XY3KB5WP	3kg/0.5g	42 071 Ft		
		XY6KA1WP	6kg/1g	42 071 Ft		
		XY7.5KA1WP	7.5kg/1g	42 071 Ft		
		XY15KA5WP	15kg/5g	42 071 Ft		
		XY30KA5WP	30kg/5g	42 071 Ft		
<b>TOUCH SYSTEM WEIGHING INDICATOR</b>						
XY-INT Indicator and Platform		XY10INTA	10kg/1g	128 429 Ft	300x400 mm	1. 7inch touch panel; 2.*Weighing,density,check weight,label print,%,counting; 3.USB input and output; 4.Information enter; 5. Connect label printer;
		XY30INTA	30kg/1g	130 643 Ft		
		XY50INTA	50kg/1g	130 643 Ft		
		XY60INTA	60kg/1g	130 643 Ft		

scale		XY10INTB	10kg/0.1g	130 643 Ft	INDICATOR	6.custom capacity and resolution;
		XY20INTB	20kg/0.1g	137 286 Ft		
		XY30INTB	30kg/0.1g	137 286 Ft		
		XY-INT indicator	CUSTOM	66 429 Ft		
<b>PRECISION WEIGHING SCALE</b>						
XY-RBG 高精度 桌称		XY6202RBG	6200g/0.01g	411 857 Ft	195x175 mm	1. EMF Sensor 2. External calibration 3. LCD with backlight with adjustable contrast 4. Plastic protection cover 5. Under hook weighing
		XY8202RBG	8200kg/0.01g	482 714 Ft		
		XY10102RBG	10100g/0.01g	522 571 Ft		
		XY12102RBG	12100g/0.01g	597 857 Ft		
<b>ECONOMY BALANCE</b>						
XY-EB Economy Balance		XY3KB1EB	3000g/0.1g	3 454 Ft	157x157mm	1. Stainless steel Pan 2. USB incharge line (no power adopter) 3. AAA dry battery x 2 4. g,oz,ct,tl 5. 24pcs/CTN 18KG G.weight/CTN
		XY6KB5EB	6000g/0.5g	3 454 Ft		
		XY10KA1EB	10000g/1g	3 454 Ft		
<b>TEXTILE BALANCE</b>						
XY-2CS/ JA-S Textile scale		XY2002CS	210g/0.01g	19 043 Ft	Φ130mm	1.Stainless steel Pan, 2.Rechargeable battery; 3.Tare, counting functions; 4. LCD wide viewing angles; 5.OZ, ct, tex, Nm, D, g/m2;OZ/yd2 7 unit conversion.
		XY3002CS	310g/0.01g	19 929 Ft		
		XY6002CS	610g/0.01g	22 143 Ft		
		XY10002CS	1000g/0.01g	24 357 Ft		
		YG101	110g/1mg	248 000 Ft	Φ90mm	* OZ, ct, tex, Nm, D, g/m2;OZ/yd2 7 kinds of textile unit conversion. Set time, date, 100 back, 200 back, 400 return sample long. Automatic storage 100 sets of data, and printable arbitrary deletion testing, each weighing values and statistical maximum, minimum.
		JA203SD	210g/1mg	248 000 Ft		
<b>DENSITY BALANCE</b>						
Density Balance		XY3002CM	310g/10mg	57 571 Ft	Pan Size: 160x70mm Tank Size: 139x100x78m m	1.Density frame 2. Electronic Balance 3.RS232 interface. 3.Solid/liquid Test *Liquid density Kits +25 USD
		XY6002CM	610g/10mg	57 571 Ft		
		JA600CM	610g/10mg	110 714 Ft		1.Density frame 2.Electronic Balance 3.RS232 interface. 3.Solid/liquid/Powder/Absorbent Test *Draft Sheild +20 USD *Liquid density Kits +25 USD (50ml)
		JA303M	310g/5mg	110 714 Ft		

		JA203M	210g/1mg	110 714 Ft		*Powder bottle + 10 USD
<b>SPARE PARTS OPTION</b>						
<b>LABEL PRINTER</b>		PRINT-S		35 429 Ft		Label Printer only for special mode option 50x60mm or 40x30mm label *label paper option
<b>PRINTER</b>		PRINT-01		22 143 Ft		Thermosensitive printer
		TXY100		93 000 Ft		Needle printer Printable header Printable time and date
		TXY110		115 143 Ft		Needle printer Printable header Printable time and date Input serial number number and product English name
<b>SPARE PARTS OPTION</b>						
<b>GRINDER</b>				11 071 Ft		Grain crusher
<b>Aluminum Box</b>				26 571 Ft		Only for XY-MW Series
<b>Temperature Meter</b>		0.1/1 °C		11 071 Ft		Only for moisture analyzer temperature calibration
<b>Sample Pan</b>		φ 90mm 50 ↑/Box		4 429 Ft		90mm only for XY-MW Series 100mm only for XY-MS Series
		φ 100mm 80 ↑/Box		8 857 Ft		
<b>Ioniser</b>		Table mode		205 929 Ft		Distance sample-ion source: About 5-40 cm Discharge time: 9s/5cm, 13s/10cm, 100s/40cm (with fan turned on) Ozone concentration: 0-0.05 ppm (2cm from ion source)
		Fix in SEM		143 929 Ft		Maximum air volume: 0.06 cm <sup>3</sup> / min Ambient conditions: 0-50°C, 20-80% air humidity (non-condensing) AC adapter (primary): AC 100-240V, 50-60 Hz Size: 11 × 10.5 × 4.5cm Weight: 310 g Installation site: device may be only used indoors

LAB-EX Laborkereskedelmi Kft.

Tel: 06-1-206-2455

E-mail: [labex@labex.hu](mailto:labex@labex.hu)

Web: [www.labex.hu](http://www.labex.hu)